User manual

MultiGlide Single Patient Use

immedia
by Etac

IMM1065_IFU - 2025-10-27 - Rev. 12.0

www.etac.com

Symbols

4

General

10

English	10
Svenska	14
Norsk	18
Dansk	22
Suomi	26
Deutsch	30
Nederlands	34
Français	38
Italiano	42
Español	46

Application suggestions

50

MultiGlide Single Patient Use	51
-------------------------------	----

Warnings	53
----------	----

Information	54
-------------	----

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Immedia MultiGlide Single Patient Use

Product label

1 PC

REF **IM302/1**

MultiGlide Single Patient Use vertical XXL

Size: W1390xL900 mm. Max. 360 kg

LOT **xxxxxxx** MD **CE**   

3

 Etac Immedia A/S - Parallelvej 3 - DK-8751 Gedved - www.etac.com

1

EN: QR-code
SV: QR-kod
NO: QR-kode
DA: QR-kode
FI: QR-koodi
DE: QR-Code
NL: QR-code
FR: Code QR
IT: Codice QR
ES: Código QR

2

EN: Back of the box: Product label
SV: Lådans baksida: Produktetikett
NO: Baksiden av esken: Produktetikett
DA: Bagside af æske: Produktmærkat
FI: Pakkauksen takasivu: Tuotteen etiketti
DE: Rückseite des Kastens: Produktetikett
NL: Achterkant van de doos: Productlabel
FR: Arrière de la boîte : Étiquette du produit
IT: Retro della confezione: Etichetta prodotto
ES: Parte posterior: Etiqueta del producto

3

EN: (01) EAN-code, (11) date, (21) Serial number
SV: (01) EAN-kod, (11) datum, (21) Serienummer
NO: (01) EAN-kode, (11) dato, (21) Serienummer
DA: (01) EAN-kode, (11) dato, (21) Serienummer
FI: (01) EAN-koodi, (11) Päivämäärä, (21) Sarjanumero
DE: (01) EAN-code, (11) Datum, (21) Serial number
NL: (01) EAN-code, (11) datum, (21) Serial number
FR: (01) Code-EAN, (11) date, (21) Numéro de série
IT: (01) Codice EAN, (11) data, (21) Numero di serie
ES: (01) Código EAN, (11) fecha, (21) Número de serie
CZ: (01) EAN kód, (11) datum, (21) Sériové číslo
SI: (01) EAN koda, (11) datum, (21) Serijska številka
PL: (01) kod EAN, (11) data, (21) numer seryjny

General

Thank you for choosing an Etac device.

This manual must be read thoroughly before you start using the device as it contains important safety information and gives instructions on how to use the device correctly. Visit www.etac.com for all relevant product documentation, for example the user manual, maintenance and recycling instructions. You can select your language via the “International” and “Local websites” link.

You can also consult the QR code on the device for further information and instructions.

Copyright

This document is the property of Etac Immedia A/S. The content shall not be modified in whole or part without prior approval of Etac Immedia A/S.

Disclaimer

Etac continuously makes improvements to our devices and information. We therefore reserve the right to modify our devices and instructions for use without prior notice. Consult www.etac.com to find the latest version of the documentation.

CE Declaration of Conformity

The device mentioned in this manual and the different models hereof comply with the Medical Device Regulation (EU) 2017/745.

The device is classified as a Class I medical device.

Product standards

The device has been tested and complies with the requirements applicable for this type of device in the below mentioned standards:

- ISO 21856:2022 Assistive products — General requirements and test methods.

- DS/EN ISO 20342-1:2019 Assistive products for tissue integrity when lying down - Part 1: General Requirements.

- The devices comply with EN 1021-1:2014 and EN ISO 12952-1:2010 resistance to ignition by cigarette (but not flame).

Warranty

Please refer to www.etac.com/support.

If the customer makes adjustments, repairs or uses combinations not predetermined by Etac, the Etac CE certification and Etac warranty will not be valid. If in doubt, please contact Etac.

Notice

Any serious incident that has occurred in relation to the device should be reported to your distributor/representative and the national competent authority within your region in a timely manner. The distributor/representative will forward the information to the manufacturer.

Pre-sales information

Please visit www.etac.com.

Device description

Etac supplies a wide selection of Positioning Horizontal devices for different types of transfer and positioning situations.

The Immedia MultiGlide Single Patient Use (SPU) is a disposable sliding sheet made from durable, low-friction nylon material. It is intended for use with a single patient to support hygienic and safe transfers and repositioning. Although designed for single-patient use, the sheet may be reused multiple times with the same user.

Using MultiGlide SPU eliminates the need for laundering and reduces storage requirements.

MultiGlide Single Patient Use are available in different sizes and models supplied in boxes containing either 10 or 20 units or on a Roll. Each sheet is individually packaged in a sealed polybag, and the compact box design allows for convenient storage and easy access within clinical departments.

See www.etac.com for a complete overview of possible variants.

Clinical benefits

The clinical benefits of the devices are:

- Reduces the need for carers or carer's time
- Reduces work-related injuries or strain of carer
- Enhance user's well-being
- Reduces risk of shear or pressure injury
- Facilitates micro-and macro movements
- Enables users to remain independent

- Facilitates assisted movement
- Relieves pain
- Reduces the risk of infection
- Enables early rehabilitation

Materials

100% nylon.

Introduction

Definitions

The person referred to as “the user” is the person who is transferred in or positioned with the device. “The carer” is the person who manages the device.

Intended purpose

The positioning horizontal are assistive devices intended for alleviation of, or compensation for, a functional impairment due to an injury or disability. The devices are designed for an user requiring reduction of friction at pressure points due to reduced mobility or physical strength.

Intended user

The intended users of the devices are based on the users' functional ability and not on a specific diagnosis, health condition or age. Ensure the user complies with the SWL/user mass and that a clinical risk assessment is conducted.

No limitation on height.

Secondary users of the devices are carers providing assistance and clinicians/technicians who set up the device.

Intended environment

The devices are intended for indoor use.

Intended application

The device is intended for 'Transient' use, meaning continuous use for less than 60 minutes.

The device is intended for use with intact skin only.

The device is intended for single patient use.

Expected service life

The expected service life of Positioning Horizontal are 2 years. The lifetime of the device may vary depending on usage frequency, loads and washing/drying frequency and conditions.

Indications

Indication for use is disability, including, but not limited to, immobile or passive users, who are in need of support to be safely and comfortably positioned and turned in lying positions or in need of mobilisation or repositioned frequently to minimise the risk of pressure ulcers and other adverse events associated with immobility. Other medical conditions are patients with pain or discomfort, patients where rolling would be dangerous, pediatrics, geriatric, patients with dementia or psychiatric problems, and obese persons.

Contraindications

There are no known contraindications.

Precautions

Do not let the device block the user's airways.

Training requirements

Only carers and/or users who have read and understood the user manual are allowed to use the device.

Safety information

The safety of both user and carer must never be jeopardized – if you suspect that the device has a malfunction you must put it in quarantine marked “out of order” and call your local distributor/representative.

Safety signs, safety measures and warnings symbols are given in the context relevant for the actual handling and/or operating situation – it is mandatory the given instruction is followed to avoid dangerous/hazardous situations.

Signs, symbols and instructions placed on this device must never be covered or removed and must remain present and clearly legible throughout the entire lifetime of the device. Immediately replace or repair illegible or damaged signs, symbols, and instructions. Contact your distributor/representative for instructions.

Patient transfer

Transfer of a user will always pose a certain risk. Therefore, a risk assessment must always be carried out before the transfer is performed.

Warnings:

- Pull/Push with a minimum of force.
- There is a risk that the user may slide off. Never leave the user alone on the edge of the bed.
- Never leave the product on the floor.
- Always remove the MultiGlide SPU after use.
- The distance between MultiGlide Single Patient Use and the edge of the bed should not be less than 15 cm (5.9 inches).

Risk assessment

Do a risk assessment of the device, the individual user and the positioning, and document your analysis. You as the carer are responsible for the safety of the user.

Plan the positioning in advance to ensure that it is performed as safe and smooth as possible.

Remember to work ergonomically according to local recommendations.

Encourage the user to assist where possible.

The device

Warning:

- Never use a faulty or damaged device as it can cause personal injury. Contact your local distributor/representative for further instructions.

Materials

All used materials which are to be in direct contact with skin are biologically assessed and evaluated according to relevant standards. However, users/carers who are allergic to any of the materials are recommend not to use the device.

Spare parts

Article numbers and spare part can be found on www.etac.com.

Combinations

When combining the device with other devices or accessories, the combination must not change the intended purpose of the devices or modify the devices in such a way that compliance with the applicable requirements may be affected. The person/company responsible for the combination of the devices must ensure that the applicable requirements are fulfilled.

Safety measures

If the device has lost low friction, dispose of it and obtain a new MultiGlide Single Patient Use for the user.

Service, cleaning and maintenance

Maintenance and service

Maintenance and inspections of the device are wholly the responsibility of the owner of the device. Failure to maintain a device in accordance with the instructions may compromise the safety of the user and/or carer and may affect the expected lifetime of the device.

Inspection of the device must be carried out regularly.

Check:

- that the materials are intact/faultless.
- never try to repair a product yourself.

Cleaning

The device cannot be washed. It can be cleaned with a solvent-free cleaning agent with a pH level between 5 and 9, or a 70% disinfection solution.

Transportation, storage and disposal

Transportation and storage

When the device is not in use, it can be stored under the following conditions:

Keep packaging away from children.

Keep away from fire.

The device can be stored and transported under temperatures between -30 to +40 °C.

Disposal

The device contains no hazardous substances that require special disposal. Check if there are options available for recycling in an environmentally friendly way and follow the national provisions that apply. Dispose of the device according to the national provisions that apply.

Recycling

See the recycling instructions on www.etac.com, and follow the national provisions that apply.

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EN: Information

01. Application suggestions, MultiGlide Single Patient Use has multiple purposes, the user's special needs and risk assessment define the final solution.
02. MultiGlide Single Patient Use is ideal be positioned under pressure points – where the user "clings" to the underlying surface.
03. The carer glides one hand between the two sheet layers, locates the corner on the opposite side and pulls it slowly towards him/ herself, turning MultiGlide Single Patient Use inside out.

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Etac Ltd.