


User manual

Hi-Loo Fixed/Supporter

etac[®]
Creating Possibilities

74949 - 2025-02-28 - Version 21

sv en da de it nl no bg cs el es et fi fr hr is lt ms pl pt sk ar


www.etac.com

MD


sv en da de it nl no bg cs el es et fi fr hr is lt ms pl pt sk ar


7 7 8 9 10 10 11 12 13 14 15 16 16 17 18 19 20 20 21 22 23 24

A


1


8030 1116-2 / 8030 1117-2 / 8030 1118-2 / *8400 1291 / *8400 1292


8030 3012-2 / *8400 1293


8030 1114-2 / 8030 1115-2


8030 3112-2 / *8400 5101


8030 1316-2 / 8030 1317-2 / 8030 1318-2 / *8400 5099 / *8400 5100


8030 3018-2 / *8400 1294

8400 1299 / 8400 1322


8400 1297 / 8400 1298

8030 1106 / 8030 1107

sv *Reservdel, endast armstöd | en *Spare part, only armrest |
 da *Reservdel, kun armlæn | de *Ersatzteil, nur Armlehne |
 it *Ricambio, solo bracciolo | nl *Reserveonderdeel, alleen armlegger |
 no *Reservdel, kun armlene | bg *Резервна част, само подлакътник |
 cs *Náhradní díl, pouze područka | el *Ανταλλακτικό, μόνο μπράτσο |
 es * Piezas de recambio, solo reposabrazos | et *Varuos, ainult
 käetugi | fi *Varaos, vain käsinoja | fr *Pièce détachée, repose-bras
 uniquement | hr *Rezervni dio, samo naslon za ruke | is *Varahlutur,
 aðeins sætisarmur | lt *Atsarginė dalis, tik rankų atrama | ms *Alat
 ganti, hanya untuk tempat letak tangan | pl *Część zamienna, tylko
 podłokietnik | pt *Peça sobressalente, apenas apoio para os braços |
 sk *Náhradný diel, iba opierka rúk | ar *قطع غيار، مسند الذراع فقط

B

1		2		3		4		5	
6		7		8		9		10	
11									

C

1

2

3

4

7

8

5

6

UDI

Bar code GS1-128

(01)XXXXXXXXXXXXXX(11)YYMMDD(21)XXXXXX

XXXXXXXXXX

XXXXXXXXXX

REF XXXXXXXXXXXX

SN XXXXXXXX

CE

MD

Etac Supply Center AB

Långgatan 12

SE-334 33 Anderstorp

YYYY-MM-DD

D






Supporter	600 mm	490 mm	365 mm	-	220/285 mm	250 mm	460 mm	4,8 kg
Supporter adj.	620/513 mm	512/336 mm	365 mm	-	220/285 mm	240/310/375 mm	460 mm	4,8 kg
Hi-Loo Fixed 6 cm	600 mm	490 mm	365 mm	210/360 mm	195/270 mm	230 mm	460 mm	4,9 kg (2,6 kg*)
Hi-Loo Fixed 10 cm	600 mm	490 mm	365 mm	210/360 mm	195/270 mm	230 mm	460 mm	5,8 kg (3,4 kg*)


sv *Reservdel, endast armstöd | **en** *Spare part, only armrest | **da** *Reservdel, kun armlæn | **de** *Ersatzteil, nur Armlehne | **it** *Ricambio, solo bracciolo | **nl** *Reserveonderdeel, alleen armlegger | **no** *Reservdel, kun armlene | **bg** *Резервна част, само подлакътник | **cs** *Náhradní díl, pouze područka | **el** *Ανταλλακτικό, μόνο μπράτσο | **es** * Piezas de recambio, solo reposabrazos | **et** *Varuosa, ainult käetugi | **fi** *Varaosa, vain käsinoja | **fr** *Pièce détachée, repose-bras uniquement | **hr** *Rezervni dio, samo naslon za ruke | **is** *Varahlutur, aðeins sætisarmur | **lt** *Atsarginė dalis, tik rankų atrama | **ms** *Alat ganti, hanya untuk tempat letak tangan | **pl** *Część zamienna, tylko podłokietnik | **pt** *Peça sobressalente, apenas apoio para os braços | **sk** *Náhradný diel, iba opierka rúk | **ar** قطع غيار، مسند الذراع فقط

	MAX kg lbs	
Supporter	150 kg	50%
Supporter adj.	130 kg	50%
Hi-Loo Fixed 6 cm	150 kg	50%
Hi-Loo Fixed 10 cm	150 kg	50%

Fe	Zn	Al	PP	PA	TPE	HDPE	

E



Hi-Loo Fixed

i


5


Supporter

i

5


6




i


8030 1106 / 8030 1107


7


F

G

H

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Thank you for choosing an Etac product. In order to avoid damage during assembly, handling and use, it is important to read this manual and save it for future reference. You can also find it at www.etac.com. You can select your language via the "International" and "Local websites" link. Here you will also find other product documentation, such as prescriber information, pre purchase guide, and reconditioning instructions.

In the manual, the user is the person sitting in the chair. The carer is the person helping the user.

Device description

Hi-Loo Fixed and Supporter are raised toilet seats, available in different heights. They are fixed mounted on the toilet and can be used with one, two or without arm supports. The arm supports are individually foldable.

Intended purpose

Hi-Loo Fixed and Supporter (hereafter also referred to as "the device" or "the product"), are medical devices intended for alleviation of or compensation for a functional impairment due to an injury or disability. The device is designed to facilitate access to the toilet for people with restricted functional ability. The raised seat and/or arm supports intend to provide stability and support when sitting down and standing up before and after performing toilet needs.

Intended user group

The user group for the device is based on the individual's functional ability and not on a specific diagnosis, health condition or age. It is intended for individuals with a height of 146 cm or more, or a mass of 40 kg or more.

Secondary users of the device are caregivers providing assistance and clinicians/technicians who set up the device.

Intended environment

The device is intended for indoor use at home environments or institutions and is suitable for use in bathrooms, but not swimming pools or similar corrosive environments.

Intended application

The device is intended for short-term and long-term use and can be applied several times per day. The device is intended to be used in contact with intact skin.

The device is intended for refurbishment and reuse

Expected service life

The expected service life is 10 years. For complete information regarding the service life of the device, see www.etac.com.

Indications

Indication for use is disability, including, but not limited to, an individual with fluctuating balance issues, physiological and/or functional impairment, and/or a history of falls.

Contraindications

There are no known contraindications.

Warnings

Warnings describing a risk element for a specific action or setting of the device can be found in the relevant section.

Ignoring the instructions for use may result in injury.

A defective device must not be used. Risk of injury.

Do not to stand on the device. Risk of injury.

Head trap between toilet water cistern and the vertical part of the arm support.

Declaration of conformity

The device meets the requirements of the Medical Device Regulation (EU) 2017/745.

The device is tested and meets the requirements of the ISO 17966.

Pre-purchase information

Pre-purchase information can be found at www.etac.com.

Adverse event

In case of an adverse event occurred in relation to use of the device, it must be reported to your local dealer and the national competent authority in a timely manner. The local dealer will forward information to the manufacturer.

Customization

is everything that goes beyond the manual's instructions and settings. Device specially adapted by the customer must not retain Etac's CE marking. Etac's warranty expires.

Device description

1. Content of the delivery
2. Seat raiser with lid
3. Arm support
4. Product label
5. Item no label

..... [Figure A](#)

Symbols [Figure B](#)

Symbols in the manual and on the device:

1. Warning, precaution or limitation.
2. Useful advice and tips.
3. Material for recycling.
4. Tools screwdriver
5. The commode must be securely anchored
6. Maximum user weight (see technical data)
- 7-8. Cleaning (see maintenance)
9. Product weight
10. Refer to user manual.
11. Refer to user manual (blue and white).

Explanation item no label: [Figure C](#)

1. Product name
2. Product description
3. Serial number
4. Item number
5. Manufacturer and date of manufacture
6. Barcode acc. to GS1-128 GTIN-14 and serial number*
7. CE marked
8. Medical device

*The device's date of manufacture can be read from the barcode on the device. The number 11 is shown below the barcode in brackets. The number combination after these brackets is the date of manufacture.

Technical data [Figure D](#)

Mounting [Figure E](#)

The toilet must be securely anchored to the floor.

For wall hung porcelain, mount armrests before step 4.

Mounting arm support [Figure F](#)

Mounting fixed arm support [Figure G](#)

Height adjustable arm supports [Figure H](#)

1. Adjusting height
2. Adjusting width
3. Folding, fixing

Maintenance

Clean the device with a solvent-free detergent which has a pH of 5-9, or with a 70% disinfectant solution. Can be decontaminated at max. 85 °C if so required.

Constituent materials are corrosion-resistant. 5-year guarantee against material and manufacturing defects. For terms and conditions, see www.etac.com.

Storage and disposal

The device should be stored indoors in a dry place at a temperature above 5 °C. If the device has been stored for a long time (more than four months), its function must be checked by an expert before use.

To be recycled according to national regulations.

Accessories

- Adjustable arm supports for Supporter
- Seat unit for angled Hi-Loo Fixed
- Cross bar
- Screw and attachment kit special Supporter
- Attachment kit

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**Etac Ltd.**
**Promefa AG**