



VELA Tango 300 & 310

NEW
Up to 300 kg
Elec. aid to get up
4XL-seat
(w:72 d:50)



VELA Tango 300

Maximum weight load: 200 kg. The user gets as standard a permanent sitting height without the possibility for adjustment.



VELA Tango 300E

Maximum weight load: 200 kg. The user can change the seat height and get up and sit down with electric height adjustment.



VELA Tango 310E

Maximum weight load: 300 kg. The user can change the seat height and get up and sit down with electric height adjustment.

Users up to 300 kg can now become more self-reliant using a robust VELA chair. This means that it will be possible to mobilize and be more active in everyday life without unnecessary fatigue. The user can push/walk around, seated, as the chair has wheels, and at the same time have a stable, safe chair with a hand-operated brake. Electric height adjustment helps the user to get up and sit down and easily change the seat height.

Obese users can mobilize themselves from a bed or chair into a safe and stable chair with a brake, which makes it easier to carry out activities. The VELA Tango 300/310 series is very robust and can be used for activities at home, at work, and during rehabilitation and training.

The user can breathe more easily and ensure oxygenation of the blood with an upright sitting position. Legs, knees, and hips are also relieved by performing the activity while sitting, which results in less muscle/joint pain and less fatigue. It is possible to push/walk the chair around while sitting and be more self-reliant in various activities, e.g. cooking, and thereby prevent immobility.

It is safe to get up and sit down, using the sturdy armrests on a secure, braked chair. There are also good conditions for comfort, as seat height, backrest, and armrests can be adjusted according to the user's physique and need for support.

Electric height adjustment can be used as an aid to get up and allows to adjust the sitting height for different activities at the touch of a button. In that way, the user can save energy during the day.

BENEFITS

- :: Robust CE marked chair for people weighing up to 300 kg
- :: Help getting up and sitting down with the electric height adjustment
- :: Be supported, stable and get a comfortable sitting position providing open airways thereby reduces fatigue
- :: Relief of the hip and knee joints and the possibility of managing one's energy by pushing/walking the chair around
- :: Feel safe and get a stable sitting position in daily activities with the central brake

STANDARD CONFIGURATION

	VELA Tango 300	VELA Tango 300E	VELA Tango 310E	
MEASUREMENTS AND WEIGHT	Height adjustment	Permanent height or 150mm movement by means of gas (optional)	150mm or 200mm movement (electrical)	150mm movement (electrical)
	Sitting height	at permanent height 42-52 cm (adjusted by means of a tool) at 150 mm gas: 44-59 cm	at 150 mm: 42-57 cm at 200 mm: 49-69 cm	at 150 mm: 42-57 cm
	Seat width x depth	58 x 48 cm (XXX seat)	58 x 48 cm (XXX seat)	72 x 50 cm (4XL seat)
	Sitting depth (adjusted by means of a tool)	46-56 cm	46-56 cm	46-56 cm
	Backrest width x height	38 x 43 cm (CX backrest)	38 x 43 cm (CX backrest)	38 x 43 cm (CX backrest)
	Height adjustment of backrest (adjusted by means of a tool)	0-15 cm	0-15 cm	0-15 cm
	Armrest padding, P.U.	7 x 35 cm	7 x 35 cm	7 x 35 cm
	Height of armrest (adjusted by means of a tool)	0-24 cm	0-24 cm	0-24 cm
	Distance between armrests (adjusted by means of a tool)	53-64 cm	53-64 cm	70-75 cm
	Frame size	54 x 54 cm	54 x 54 cm	54 x 54 cm
	Chair's weight	29 kg	37 kg	42 kg
	Maximum weight load	200 kg (gas 175 kg)	200 kg	300 kg
	Item number	5024	5105	5106

All VELA chairs are CE marked and comply with current standards.

VELA bariatric chairs are supplied as standard in black Phoenix fabric, which is durable, breathable and comfortable.

We offer a large selection of artificial leather or the option of different seat fabrics if the materials need to be even more cleaning friendly.

ACCESSORIES AND CUSTOMISED SOLUTIONS

VELA Tango 300/310 is exceptionally durable and can withstand constant use under heavy loads.

The chair can be used by multiple users in rehabilitation and day centres

Complete accessory list with item numbers, pictures and descriptions can be found at www.vela.eu.

WORK ENVIRONMENT

- :: The helpers work environment is improve with electric height adjustment, as it supports the user from sitting to standing position.
- :: Easy to help the user level with a bed or chair, e.g. for transfers, due to electric height adjustment.
- :: Possible to order a push handle if the user needs to be moved over short distances.



VELA Tango 300/310 features