

ErgoMulti Table - Height Adjustable Training Table

Product details

↑ Manual height adjustment	34 cm (56-90 cm)
🔧 Operation	Handle
↗ Dimensions	H: 56-90 cm x W: 90 cm x D: 70 cm (Inner width: 73.6 cm)



Description

The ErgoMulti Table is a height adjustable training table ideal for therapy sessions where two people are sitting facing each other. The two ergonomic curves on the table sides provide for a healthy sitting position for a user and a therapist at the same time.

As there are no cross bars under the table top, wheelchair users have no problem getting close to the table and into a comfortable position. The training table is height adjustable from 56-90 cm by means of a detachable handle.

The standard version of the table comprises a boomerang system which allows easy change of the cut out sizes between 30/40 cm diameter or using the half moon completely fill the cut out

Specifications

Item number	20-70410
Power	Manual
Operation	Handle
Configuration	Two chamfered corners Incl. 2 boomerangs and 1 half moon
Max. Load on the table (static)	100 kg
Max. Load during lifting (dynamic)	50 kg

Item No: 20-70410

ErgoMulti Table - Height Adjustable Training Table

Materials Steel St. 37. Various plastic components. Tabletop: Laminate

Surface treatment Powder coating. Chromate

Standard colour Light grey

Included Brake wheels

CE certified

Line drawing

