



VELA Meywalk 4



VELA Meywalk 4 Large



VELA Meywalk 4 Medium



VELA Meywalk 4 Small



VELA Meywalk 4 Micro

VELA Meywalk 4 is a series of gait trainers designed to encourage self-movement in children and adults suffering from motor disabilities. This easy-to-use gait trainer provides optimal support to help the user stand and walk independently. The ability to move freely and interact socially at eye level promotes mobility and social participation, enhancing the user's mental and physical wellbeing.

VELA Meywalk 4 promotes self-movement while stimulating and supporting the user's natural pattern of movement. Designed with a dynamic seat and a unique spring system, the gait trainer follows the user's movements, ensuring optimal weight-bearing. The body ring allows for movement of the trunk while supporting an optimal upright posture.

VELA Meywalk 4 is safe and easy to access. Remember to apply the brake on the gait trainer and use the integrated swing bar to lower the seat, and then the user can relocate to the gait trainer with minimal assistance. The seat is easily raised using the built-in tilt bracket.

The user is not dependent on their forearms for weight-bearing, promoting social participation and enabling the person to eat/drink while using the gait trainer. The integrated hand bar is easily adjusted to fit the user if additional support is needed.

Daily standing in an upright position influences joint development, bone density and muscle strength profoundly. So, mobility is maintained, while endurance and balance increase. Developing the ability to stand and walk with a VELA Meywalk 4 vastly improves blood circulation, digestion and respiration, promoting natural appetite, better sleep and overall well-being, which is vital and often a prerequisite for managing others daily activities.

BENEFITS

- :: VELA Meywalk 4 stimulates and supports the body's natural pattern of movement.
- :: User can be standing/walking despite the lack of or reduced standing function
- :: Grows with the user thanks to the variable height adjustment of seat and body ring.
- :: Versatile and adjusts easily without the use of tools.
- :: Easy to separate for transport purposes.

STANDARD CONFIGURATION

MEASUREMENTS AND WEIGHT		VELA Meywalk 4 Micro	VELA Meywalk 4 Small	VELA Meywalk 4 Medium	VELA Meywalk 4 Large
	Length	77 cm	85 cm	86 cm	100 cm
	Width, outer measures	51 cm	62 cm	62 cm	68 cm
	Width, inner measures	42 cm	49 cm	50 cm	57 cm
	Maximum turning width	95 cm	105 cm	106 cm	121 cm
	Seat width:				
	Front	7 cm	9 cm	10 cm	10 cm
	Middle	5,2 cm	7,5 cm	7,2 cm	7,2 cm
	Back	11,5 cm	13,5 cm	15,5 cm	15,5 cm
	Seat depth	24 cm	27 cm	32,5 cm	32,5 cm
	Seat height/stride length*	25-45 cm	44-67 cm	60-84 cm	80-104 cm
	Height of body ring*	52-65 cm	74-99 cm	101-133 cm	120-152 cm
	Size of body ring (select size based on waist measurement of the user)	70 or 90 cm	70, 90 , 105 or 120 cm	70, 90, 105 , 120 or 130 cm	70, 90, 105, 120 or 130 cm
	Weight (without accessories)	13 kg	25 kg	32 kg	34 kg
	Maximum weight load	60 kg	100 kg	130 kg	130 kg
	Weight interval (colour on springs) (Select based on the body weight of the user)	10-15 kg (orange) 15-20 kg (purple) 20-50 kg (white) 45-65 kg (yellow)	20-50 kg (white) 45-65 kg (yellow) 65-80 kg (blue)	45-65 kg (yellow) 65-80 kg (blue) 80-100 kg (red) 100-130 kg (silver)	45-65 kg (yellow) 65-80 kg (blue) 80-100 kg (red) 100-130 kg (silver)
	Colour of the gait trainer	Black	Green	Red	Black
	Frame	Aluminium	Steel	Steel	Steel
	Expected lifetime	5 years	5 years	5 years	5 years
	Brake	Yes	Yes	Yes	Yes
	Walking resistance brake	Yes	Yes	Yes	Yes
	Backstop	Yes	Yes	Yes	Yes
	Item number	MS014000	MS013000	MS011200	MS011250

* Height is measured with half compression of springs.

Standard is in **bold**

The gait trainer is not intended to be folded or dismantled.

The user manual can also be found on vela.eu homepage. The user manual must follow the gait trainer. If the user is not able to read or cognitive able to understand the user manual a helper must step in. The gait trainer may only be used under supervision.

All VELA Meywalk 4 gait trainers are CE marked and comply with current standards.

The gait trainer is approved for indoor use.

The standard gait trainer is supplied with upholstered parts in artificial leather and polyurethane foam which are durable and easy to clean. The gait trainer is not intended to be disinfected.

ACCESSORIES AND CUSTOMISED SOLUTIONS

All VELA Meywalk 4 gait trainers are versatile and customisable to the user's specific needs thanks to a wide range of accessories, for example, a hip and trunk support to optimise the body posture.

Distributor Details:

You can find a complete accessory list with item numbers, pictures and descriptions at vela.eu

Do only use original accessories from VELA.

WORK ENVIRONMENT

The VELA Meywalk 4 gait trainer is versatile and adjusts easily without tools, especially beneficial in healthcare environments where several users may share one gait trainer. The gait trainer is not intended to be used in combination with other products – due to dangerous combinations and non-flame-resistant materials.

The seat can be lowered using the integrated swing bar so users can relocate to the gait trainer with minimal assistance. It is also possible to lift the user into the Meywalk 4 gait trainer.